

## Product Manual

<b>Product Name</b>	Pfu DNA Polymerase
<b>Source</b>	Recombinantly expressed in Escherichia coli
<b>Catalog Number</b>	CSB-DEM035
<b>Physical Form</b>	Liquid
<b>Enzyme Activity</b>	1 U/ $\mu$ L
<b>Storage Conditions</b>	-20 $\pm$ 5 $^{\circ}$ C
<b>Storage Buffer</b>	10 mM Tris-HCl, 100 mM KCl, 1 mM DTT, 0.1 mM EDTA, 0.5% Tween 20, 0.5% NP-40, 50% Glycerol
<b>Activity Definition</b>	The amount of enzyme required to incorporate 10 nmol of deoxynucleotides into acid-insoluble material in 30 minutes at 72 $^{\circ}$ C is defined as 1 unit (U).
<b>Quality Control</b>	Free of exonuclease and endonuclease contamination
<b>Shelf Life</b>	24 months

### Product Description

Pfu DNA Polymerase is a highly heat-stable DNA polymerase derived from the thermophilic bacterium Pyrococcus furiosus and recombinantly expressed in Escherichia coli. It exhibits high amplification efficiency and fidelity. Due to its strong 3' - 5' exonuclease activity, it possesses higher fidelity compared to Taq DNA Polymerase, along with efficient amplification capabilities at a speed of up to 1 kb/10s. Pfu DNA Polymerase generates blunt-ended products in the amplification reaction.

### Product components

Component NO.	Component Name	Specifications		
		0.6mL	3mL	6mL
1	5 $\times$ Pfu buffer	0.6mL	3mL	6mL
2	Pfu DNA Polymerase	100 U	500U	1000 U

WUHAN HUAMEI BIOTECH CO.,LTD

📍 No.818 Gaoxin Avenue, Wuhan Hi-tech Medical Devices Park, Donghu High-tech Development Zone 430206, Wuhan City, Hubei Province, P.R. China.

🌐 [www.cusag.cn](http://www.cusag.cn) / [www.cusagivd.com](http://www.cusagivd.com)

✉ [cusag@cusag.cn](mailto:cusag@cusag.cn)

☎ +86-27-65521556/+86-27-87196282 Ext 853

📞 +86-27-87196150

## Usage of Product Components

### Recommended Reaction System

Components	Addition amount
ddH <sub>2</sub> O	Up to 30 $\mu$ L
5 $\times$ Pfu buffer	6 $\mu$ L
10 mM dNTPs	0.6 $\mu$ L
Forward Primer (10 $\mu$ M)	0.4-1 $\mu$ L
Reverse Primer (10 $\mu$ M)	0.4-1 $\mu$ L
Pfu DNA Polymerase	1 $\mu$ L
Template DNA	X $\mu$ L

#### a. Primer Concentration:

Generally, a primer concentration of 0.2  $\mu$ M in the reaction system yields good results.

Adjust the primer concentration within the range of 0.1  $\mu$ M to 1.0  $\mu$ M if the reaction performance is poor.

#### b. Template Concentration:

Genomic DNA from animals and plants: 0.1-1  $\mu$ g, E. coli genomic DNA: 10-100 ng,  $\lambda$ DNA: 0.1-10 ng, plasmid DNA: 0.1-10 ng. If using undiluted cDNA, the volume should not exceed 1/10 of the total qPCR reaction volume.

#### c. Polymerase Concentration:

WUHAN HUAMEI BIOTECH CO.,LTD

📍 No.818 Gaoxin Avenue, Wuhan Hi-tech Medical Devices Park, Donghu High-tech Development Zone 430206, Wuhan City, Hubei Province, P.R. China.

🌐 [www.cusag.cn](http://www.cusag.cn) / [www.cusagivd.com](http://www.cusagivd.com)

✉ [cusag@cusag.cn](mailto:cusag@cusag.cn)

☎ +86-27-65521556/+86-27-87196282 Ext 853

☎ +86-27-87196150

The enzyme amount can be adjusted between 0.25-1  $\mu$ L. Increasing the enzyme amount generally improves the amplification yield but may decrease specificity.

**d. Mg<sup>2+</sup> and Additives:**

The Mg<sup>2+</sup> concentration in most PCR reaction systems should be within the range of 1.0-5.0 mM. However, for difficult amplification samples such as high GC content DNA, additives like DMSO or formamide may need to be added to the PCR reaction system.

**Recommended PCR Reaction Program**

Temperature	Time	Number of cycles
95°C	3 min	1
95°C	10 s	25-35
(T <sub>m</sub> -5)°C	20 s	
72°C	1 kb/min	
72°C	5 min	1

**e. Annealing Temperature and Time:**

The annealing temperature should be adjusted based on the primer's T<sub>m</sub> value, generally set 3-5°C below the primer's T<sub>m</sub> value. The recommended annealing time is 20 seconds, adjustable within the range of 10-30 seconds