

Product Manual

| | |
|---------------------|--|
| Product Name | 2×Multiplex ARMS qPCR Mix |
| Catalog Number | CSB-DKT075 |
| Physical Form | Liquid |
| Storage Conditions | -20 ±5°C, avoid repeated freeze-thaw cycles; can be stored at 2~8°C in the dark for frequent use. |
| Shipping Condition | ≤0°C; dry ice transportation |
| Product Composition | Buffer system, dNTPs, ARMS Taq DNA Polymerase |
| Quality Control | All components have been tested and found to be free of nucleases, endonucleases, and RNase residues |
| Shelf Life | 12 months |

Product Description

2× Multiplex ARMS qPCR Mix is a reagent kit developed based on the principle of amplification blocking. It has strong blocking ability and can effectively perform gene typing mediated by ARMS-PCR. The buffer and Taq enzyme are specifically optimized and contain the necessary ions and dNTPs for amplification, eliminating the need for additional additives. This system can detect various gene mutations, detecting mutant alleles as low as 0.05% in genomic DNA and 0.01% in plasmid DNA.

Product Components

| Component No. | Component Name | Specification | Specification | Specification |
|---------------|---------------------------|---------------|---------------|---------------|
| | | (50T) | (100T) | (500T) |
| 1 | 2×Multiplex ARMS qPCR Mix | 0.5mL | 1mL | 5mL |

Operating Instructions

Recommended Reaction System

| Components | Additions |
|--------------------------------|------------------|
| ddH ₂ O | Up to 30 μ L |
| 2×Multiplex ARMS qPCR Mix | 15 μ L |
| Upstream Primer(10 μ M) | 1 μ L |
| Downstream Primer(10 μ M) | 1 μ L |
| Fluorescence Probe(10 μ M) | 0.3 μ L |
| TemplateDNA | X μ L |

[Note]:

a. Primer Concentration: Generally, a primer final concentration of 0.2 μ M in the reaction system yields good results. If the reaction performance is poor, the primer concentration can be adjusted within the range of 0.1 μ M to 1.0 μ M.

b. Template Concentration: Genomic DNA from animals and plants: 0.1-1 μ g; Escherichia coli genomic DNA: 10-100 ng; λ DNA: 0.1-10 ng; Plasmid DNA: 0.1-10 ng. If the template is undiluted cDNA, the volume used should not exceed 1/10 of the total qPCR reaction volume.

Recommended PCR Reaction Program

| Temperature | Time | Cycles |
|-------------|---------|--------|
| 95°C | 3 min | 1 |
| 95°C | 10 sec | 40-45 |
| 48-68°C | 30 sec* | |

[Note]:

c. Annealing Temperature and Time: The annealing temperature should be adjusted based on the primer's T_m value, generally set to 3-5°C below the primer's T_m value.

* Collect fluorescence.