

## Product Manual

<b>Product Name</b>	HiFi PCR Kit
<b>Catalog Number</b>	CSB-DKT039
<b>Formation</b>	Liquid
<b>Storage</b>	-20 ±5 °C, avoid repeated freeze-thaw cycles; for frequent use, store at 2~8°C, protected from light
<b>Transportation</b>	0°C or less; Transport on dry ice
<b>Contents</b>	buffer、dNTPs、HiFi DNA Polymerase、GC Enhancer
<b>Quality Control</b>	All components have been tested and confirmed to be free of nucleases, both exonucleases and endonucleases, as well as RNase residues.
<b>Validity Period</b>	12 months

### Product Description

HiFi PCR Kit is a ready-to-use premixed solution. The 2\*HiFi PCR Mix contains HiFi DNA Polymerase, dNTPs, and an optimized buffer system.

HiFi DNA Polymerase is a highly heat-stable DNA polymerase, obtained through recombinant expression in Escherichia coli. It is a next-generation high-fidelity DNA polymerase with high amplification efficiency and fidelity. Its strong 3' - 5' exonuclease activity provides higher fidelity than Taq DNA Polymerase, along with efficient amplification. When combined with the optimized system, it can achieve a speed of 1 kb/10 s and has excellent long fragment amplification ability.

### Product Components

Label	Components	Specification (50T)	Specification (100T)	Specification(500T)
1	2*HiFi PCR Mix	0.75mL	1.5mL	5*1.5mL
2	5*GC Enhancerr	0.6mL	1.2mL	5*1.2mL

### Operating instructions

#### Recommended Reaction System

Component	Volume
ddH2O	Up to 30 µL
2*HiFi PCR Mix	15 µL
Forward Primer (10µM)	1 µL
Reverse Primer (10µM)	1 µL

Template DNA 10ng-1 $\mu$ g

**【Notes】 a. Primers:** The final concentration of primers in the reaction system is 0.2 $\mu$ M to obtain better results. When the reaction performance is poor, the primer concentration can be adjusted within the final concentration range of 0.1 $\mu$ M-1.0 $\mu$ M.

**b. Template :** Animal and plant genomic DNA 0.1 ~ 1 $\mu$ g; Escherichia coli genomic DNA 10-100 ng;  $\lambda$ DNA 0.1 ~ 10ng; Plasmid DNA 0.1-10 ng. If the template is an undiluted cDNA stock solution, the volume used should not exceed 1/10 of the total volume of the qPCR reaction.

**c. GC Enhancer:** If the template is complex or rich in GC, GC Enhancer can be added to the reaction system. The storage concentration of GC Enhancer is 5 $\times$ , and the working concentration can be adjusted between 0.5 $\times$  and 2 $\times$ .

**Recommended PCR Reaction Procedure**

Temperature	Time	Cycle Number
95 $^{\circ}$ C	3 min	1
95 $^{\circ}$ C	10 sec	25-35
(T <sub>m</sub> -5) $^{\circ}$ C	20 sec	
72 $^{\circ}$ C	15-30 sec/kb	
72 $^{\circ}$ C	5 -10 min	1

**【Notes】 d. Annealing temperature and time:** The annealing temperature needs to be adjusted according to the T<sub>m</sub> value of the primer, and is generally set to be 3 to 5 $^{\circ}$ C lower than the T<sub>m</sub> value of the primer. The recommended annealing time is set to 20 sec, which can be adjusted within 10-30 sec.