

## Product Manual

<b>Product Name</b>	T4 DNA Ligase
<b>Source</b>	Recombinant expression in Escherichia coli
<b>Catalog Number</b>	CSB-DEM049
<b>Physical Form</b>	Liquid
<b>Storage Conditions</b>	-20 ± 5°C
<b>Molecular Weight</b>	55.3 kDa
<b>Quality Control</b>	No residual exonucleases or endonucleases.
<b>Shelf Life</b>	24 months

### Product Description

T4 DNA Ligase is an enzyme that catalyzes the formation of phosphodiester bonds between the 5'-P end of one DNA or RNA molecule and the 3'-OH end of another DNA or RNA molecule, either blunt-ended or sticky-ended. This catalytic reaction requires ATP as a cofactor. Additionally, T4 DNA Ligase can repair single-stranded nicks on double-stranded DNA, double-stranded RNA, or DNA/RNA hybrids. This product is suitable for connecting restriction enzyme-digested fragments, linkers, or adapters, as well as for nick repair and Ligase-mediated RNA detection.

### Product components

Component No.	Component Name	Specifications		
1	T4 DNA ligase (400U/μL)	40KU	80KU	400KU
2	10× T4 DNA Ligase Buffer	300μL	600μL	1.5mL*2

## Operating instructions

### 1. Prepare the following reaction system in a sterile microcentrifuge tube

Composition	Usage
ddH <sub>2</sub> O	Up to 20 $\mu$ L
10 $\times$ T4 DNA Ligase Buffer	2 $\mu$ L
T4 DNA ligase (400U/ $\mu$ L)	1 $\mu$ L
Vector DNA	Approximately 50-100 ng
Insert fragment	Approximately 3 times the molar amount of the vector

#### Note:

When connecting the blunt-ended vector with DNA fragments, the vector should be dephosphorylated first to prevent self-ligation. To improve ligation efficiency, 2  $\mu$ L of 50% PEG 4000 can be added to every 20  $\mu$ L reaction system.

### 2. Perform the reaction overnight at 16°C.

### 3. Transformation Experiment:

- (1) Add the ligation product to 100  $\mu$ L of competent cells (the amount of ligation product should not exceed 1/10 of the competent cells) Gently mix by flicking and incubate on ice for 30 minutes.
- (2) Heat shock at 42°C for 90 seconds, then immediately place on ice for 2-3 minutes.
- (3) Add 900  $\mu$ L of LB or SOC medium to the centrifuge tube and shake at 37°C and 220 rpm for 1 hour.

- (4) Centrifuge at 2500g for 5 minutes, remove 900  $\mu$ L of supernatant, resuspend the pellet in the remaining medium, and evenly spread the remaining cells on an appropriate antibiotic-resistant agar plate using a sterile spreader Invert the plate and incubate at 37°C overnight.